

METHODOLOGY

Open Access



Optimization of immune receptor-related hypersensitive cell death response assay using agrobacterium-mediated transient expression in tobacco plants

Sung Un Huh*

Abstract

Background: The study of the regulatory mechanisms of evolutionarily conserved Nucleotide-binding leucine-rich repeat (NLR) resistance (R) proteins in animals and plants is of increasing importance due to understanding basic immunity and the value of various crop engineering applications of NLR immune receptors. The importance of temperature is also emerging when applying NLR to crops responding to global climate change. In particular, studies of pathogen effector recognition and autoimmune activity of NLRs in plants can quickly and easily determine their function in tobacco using agro-mediated transient assay. However, there are conditions that should not be overlooked in these cell death-related assays in tobacco.

Results: Environmental conditions play an important role in the immune response of plants. The system used in this study was to establish conditions for optimal hypersensitive response (HR) cell death analysis by using the paired NLR RPS4/RRS1 autoimmune and AvrRps4 effector recognition system. The most suitable greenhouse temperature for growing plants was fixed at 22 °C. In this study, RPS4/RRS1-mediated autoimmune activity, RPS4 TIR domain-dependent cell death, and RPS4/RRS1-mediated HR cell death upon AvrRps4 perception significantly inhibited under conditions of 65% humidity. The HR is strongly activated when the humidity is below 10%. Besides, the leaf position of tobacco is important for HR cell death. Position #4 of the leaf from the top in 4–5 weeks old tobacco plants showed the most effective HR cell death.

Conclusions: As whole genome sequencing (WGS) or resistance gene enrichment sequencing (RenSeq) of various crops continues, different types of NLRs and their functions will be studied. At this time, if we optimize the conditions for evaluating NLR-mediated HR cell death, it will help to more accurately identify the function of NLRs. In addition, it will be possible to contribute to crop development in response to global climate change through NLR engineering.

Keywords: NLR-type immune receptor, Pathogen effector, Cell death, Paired NLR, Tobacco, Agrobacterium, Hypersensitive response

Introduction

Recent climate change is changing the map of crop production. An increase in temperature can aggravate crop damage caused by pathogens [1, 2]. To actively respond to these rapid changes, a rapid response through genetic engineering is required.

*Correspondence: sungun@kunsan.ac.kr

Department of Biological Science, Kunsan National University, Gunsan 54150, Republic of Korea



Plant pathogens inject effector proteins by manipulating cellular processes via type three secreted systems (T3SEs) into host cells [3]. Nucleotide-binding leucine-rich repeat receptor (NLR) proteins recognize specific effectors and trigger effector-triggered immunity (ETI), which is often associated with localized programmed cell death, known as the hypersensitive response (HR) and limits pathogen proliferation [4]. Plant NLRs are structurally and evolutionarily similar to animal NLRs [5]. How plant NLRs signal immune responses remains largely unknown when compared to similar animal systems.

Like animals, paired NLRs exist in plants, and their functions are divided into sensor NLR and helper NLR. In particular, sensor NLRs have a role in recognizing effectors via the integrated domain and inhibit the autoimmune activity of helper NLRs [6]. For example, this is the case with RPS4/RRS1. RRS1 is known to recognize effectors secreted by three different pathogens [7, 8]. Interestingly, artificially extending the N- or C-terminus of RRS1 activates RPS4 autoimmunity without effector recognition [9]. That is, an autoimmune system and an effector recognition system using RPS4/RRS1 paired NLR can be studied together.

Agrobacterium-mediated transient analysis is widely used to evaluate the function of NLRs. This is because the clear phenotype of HR cell death can be assessed with various combinations of NLRs and effectors in a relatively short time without obtaining the transgenic plants. It is known that NLR-mediated cell death is influenced by environmental factors such as temperature and humidity [10, 11]. In particular, the immunity of plants decreases with increasing temperature conditions and affects effector recognition of NLRs [12]. If so, it is important to clarify the conditions for an effective HR cell death assay using an *agrobacterium*-mediated transient system in tobacco plants.

In this study, we confirmed that humidity and plant leaf position, excluding temperature conditions, were important for effectively confirming HR cell death assays in tobacco plants. The RPS4/RRS1-mediated cell death was activated by AvrRps4 perception in tobacco plants; this HR cell death was suppressed at 65% humidity. However, at 10% humidity, HR cell death was enhanced. Additionally, the infiltrated-leaf position of tobacco is important. Position #4 of recent leaves from the tops in 4–5 weeks old plants showed the most effective HR cell death. Under these distinct conditions, performing NLR-mediated HR upon effector recognition or autoimmune activation may yield more accurate results.


Results

The NLR RPS4/RRS1 pair can be used in the cell death assay system, which is divided into autoimmune activity and effector-triggered hypersensitive cell death response

Two NLR immune receptor proteins encoded by genetically linked genes function together as paired NLRs. They are divided into sensor NLR and helper NLR according to their function. Sensor NLR usually has a non-canonical domain called the integrated domain (ID) which can recognize cognate effector proteins [13, 14]. Overexpression of helper NLR caused HR-like cell death in the absence of sensor NLR. If coexpressed sensor/helper NLRs, the HR cell death is inhibited (Fig. 1A) [13, 15]. When the N-terminus or C-terminus of sensor NLR RRS1 was extended, RPS4-dependent autoimmunity appeared without effector perception (Fig. 1A). This is probably because the protein conformational change of RRS1 loses regulation of RPS4 autoimmune activity [9]. In the engineering of sensor NLRs, it must be careful when engineering to employ specific proteins. TIR (Toll-like, Interleukin-1 receptor, Resistance protein) or CC (coiled-coil) of the helper NLR act as a signaling domain. The TIR^{RPS4} domain alone can activate TIR-dependent autoimmune responses and TIR^{RRS1} can suppress this HR cell death [16]. TIR domain mutant RPS4(SH/AA) fail to activate HR, suggesting that TIR^{RPS4} oligomerization is required for generating the HR signal (Fig. 1B). RRS1-R from accessions Ws-2 forms an immune receptor complex with RPS4 that recognizes *Pseudomonas syringae* AvrRps4, *Ralstonia solanacearum* PopP2, and unknown effectors [7, 8, 17, 18]. Thus, the paired RPS4/RRS1-R can prevent infection by three distinct pathogens and it can be used as a very useful tool for crop engineering (Fig. 1C). We summarize the following experimental procedures based on the RPS4/RRS1 system capable of performing these three types of HR cell death assays using the *Agrobacterium*-mediated transient expression in tobacco plant (Fig. 1D).

Both RRS1-S and RRS1-R fused with fluorescent proteins under low humidity conditions exhibited RPS4-dependent HR cell death

Sensor NLR is a good material for developing crops that recognize various pathogens by introducing new IDs. The positioning may alter the overall sensor NLR protein size or fuse new IDs into the N- or C-terminus [19]. Bimolecular fluorescence complementation (BiFC) methods to confirm intra-/intercellular interactions should evaluate their function when applied to sensor NLRs [9, 20]. We found that overexpression of *cCFP-RRS1-S-nCerulean* activated RPS4-dependent HR cell death in low humidity conditions, but this cell death was suppressed in the high humidity condition (Fig. 2A). At 4 dpi, we could not



detect any difference between *cCFP-RRS1-S-nCerulean/RPS4-HA* mediated autoimmune activity and *cCFP-RRS1-S-nCerulean/RPS4-HA* mediated HR cell death upon AvrRps4 perception (Fig. 2A). As a result, it was confirmed that the humidity condition at a temperature of 22 °C affected HR cell death.

In the case of RRS1-R, unlike RRS1-S, 83 amino acids are extended at the C-terminus, which is known to play a decisive role in recognizing PopP2 [21]. We tested changes in autoimmunity under different humidity conditions using RRS1R-cCFP and RRS1-R-nCerulean. Consistent with *cCFP-RRS1-S-nCerulean*, both *RRS1-R-cCFP* and *RRS1-R-nCerulean* activated RPS4-dependent HR cell death under low humidity. High humidity can inhibit HR cell death in *N. tabacum* (Fig. 2B). Thus, humidity is an important determinant of HR cell death activation by paired NLR autoimmune.

Agro-infiltrated leaf location in the tobacco plant is also a critical factor in HR cell death assays

HR cell death assays using tobacco plants often show differences in the intensities of HR cell death. In particular,

if agro-infiltration is performed using various leaves in one plant, the same cell death cannot be obtained. To confirm that these differences occurred, leaf positions of 4–5 week old *N. benthamiana* were numbered from top to bottom. In leaf position #4, coexpression of *RRS1-R-cCFP/RPS4-Myc* or *RRS1-R-nCerulean/RPS4-Myc* exhibited strong HR cell death but not leaf position #5 (Fig. 3A). Consistently, this autoimmune activity was suppressed by high humidity in *N. benthamiana* (Fig. 3B). These results suggested that HR cell death is also affected by infiltrated leaf position and can be more easily detected under a relatively low humidity condition.

Agro-infiltrating leaf position is an important factor to evaluate HR cell death in both TIR^{RPS4}-mediated cell death and RPS4/RRS1-mediated effector-triggered cell death

NLRs are divided into two groups depending on the type of N-terminal domain. The N-terminal coiled-coil (CC) domain is called CNLs (CC-NLRs), and those with the N-terminal Toll/interleukin-1 receptor (TIR) domain is called TNLs (TIR-NLRs) [22]. Evolutionally, bacterial





Fig. 2 RPS4/RRS1 autoimmune activity is enhanced by low humidity. **A** *cCFP-RRS1-S-nCerulean* triggered RPS4-dependent autoimmune cell death in *N. tabacum*. *35S::cCFP-RRS1-S-nCerulean* or *35S::cCFP-RRS1-S-nCerulean/35S::RPS4-HA* or *35S::cCFP-RRS1-S-nCerulean/35S::RPS4-HA/35S::AvrRps4-GFP* or *35S::cCFP-RRS1-S-nCerulean/35S::RPS4-HA/35S::AvrRps4^{E187A}-GFP* were transiently co-expressed in *N. tabacum* leaves. HR cell death was confirmed on the 4 days. **B** *RRS1-R-nCerulean* and *RRS1-R-nCFP* triggered RPS4-dependent autoimmune cell death in *N. tabacum*. *35S::RRS1-R-nCerulean* or *35S::RRS1-R-nCerulean/35S::RPS4-Myc* or *35S::RRS1-R-cCFP* or *35S::RRS1-R-cCFP/35S::RPS4-Myc* or *35S::RRS1-R-HF* or *35S::RRS1-R-HF/35S::RPS4-Myc* were transiently co-expressed in *N. tabacum* leaves. HR cell death was confirmed on the 4 days

TIR domain proteins have NADase enzymatic activity that generates a non-canonical variant cyclic ADPR (cADPR) molecule and cleaves NAD⁺ (nicotinamide adenine dinucleotide) [23]. The TIR domain of plant NLR

also has NADase activity and is required downstream signaling [24, 25]. The overexpression of the TIR^{RPS4} domain is sufficient to activate autoimmune responses, and oligomerization of the TIR^{RPS4} domain is required




for plant immune signaling [16]. In TIR-mediated immunity studies, the importance of TIR function evaluation is increasing.

N. benthamiana is an excellent system for studying cell death induced by the TIR/CC domains. TIR^{RPS4} was used to find suitable leaf positions for HR cell death assay in *N. benthamiana*. When co-delivered *35S::TIR^{RPS4}-GFP/35S::GFP* control in leaf positions #4-#7, the HR cell death is detected in leaf position #4 and #5 (Fig. 4A). In leaf position #4, TIR^{RPS4}-mediated autoimmune cell death was more potent than in leaf position #5 (Fig. 4A). In leaf positions #6 and #7, HR cell death was not detected (Fig. 4A). If co-expressed *35S::TIR^{RPS4}(SH/AA)-GFP/35S::GFP* control, It cannot

induce HR cell death because of inhibition of TIR^{RPS4} oligomerization (Fig. 4A).

It is also checked whether RPS4/RRS1 is functional in *N. benthamiana* system. When co-delivered *35S::RRS1-R-HF/35S::RPS4-Myc/35S::AvrRps4-GFP*, clear HR cell death was observed in the leaf positions #4, but weak HR cell death was exhibited in the leaf positions #5 (Fig. 4B). However, RPS4/RRS1-mediated HR cell death upon AvrRps4 perception is not observed in leaf positions #6-#7 (Fig. 4B). To determine whether HR intensity was determined by leaf position in *N. tabacum* as in *N. benthamiana*. Similarly, a combination of *35S::RRS1-R-HF/35S::RPS4-Myc/35S::AvrRps4-GFP* was co-expressed at leaf positions #4-#5 in *N. tabacum*. HR




cell death phenotype was checked at 4 dpi. As expected, HR cell death was stronger in leaf position #4 than #5 (Fig. 4C). In effector-triggered HR cell death, the position of the leaf used for agro-infiltration affects the intensity of HR cell death, in both *N. benthamiana* and *N. tabacum*.

Based on the results, we have summarized the important factors in the tobacco-based agro-infiltration cell death assay (Fig. 5). In most reported cases, high temperature suppressed HR as well as disease resistance. The high temperature suppressed plant resistance under various experimental conditions [10, 11]. In the HR cell death analysis, humidity and leaf position used for infiltration was the most important factors. Ultimately, it is expected to be of great help in experiments to confirm the cell death phenotype that exhibits various NLRs under these optimal conditions. In addition, when preparing a sample for co-immunoprecipitation (co-IP) or western blot analysis rather than a cell death assay, if the humidity is high

and the sample is collected from a location other than location #4, sufficient experiments are possible.

Discussion

In this study, we developed and optimized an agro-mediated transient cell death assay in tobacco plants. Tobacco plants are the most effective functional research system as various types of plant NLRs have been discovered through whole genome sequencing (WGS) and resistance gene enrichment sequencing (RenSeq) [26–31]. In particular, it is possible to quickly and easily identify which effector NLR recognizes or which domain regulates NLR autoimmune activity through phenotypes such as HR cell death. On the other hand, it is difficult to confirm the intracellular localization of the NLR fused with a fluorescent protein due to autoimmune activity. Some NLRs do not produce distinct HR responses in tobacco plants. It will be possible to determine the optimal HR cell death phenotype in NLR studies using low humidity conditions and infiltrated leaf positions.



The regulation mechanism of NLR-mediated humidity-sensitive HR is not fully understood. However, this is probably similar to temperature-sensitive HR cell death. Although there are exceptions and not much research has been done, high temperatures and humidity suppressed HR and affected plant immunity [10, 11, 32, 33]. For example, dwarf phenotype of autoimmune *snc1-1* is suppressed at 28 °C [34]. *Arabidopsis* U-box ubiquitin ligase SAUL1 regulated senescence and cell death. The *saul1-1* mutant showed autoimmune activity. An autoimmune phenotype of *saul1-1* was rescued by higher relative humidity and higher temperatures. The *saul1-1* phenotype regulated by EDS1/PAD4 dependent signaling pathway [32]. This implied that NLR-mediated cell death also might be connected with EDS1/PAD4 pathway. Salicylic acid (SA) may play an important role in the thermoregulation of plant NLR-mediated cell death [33].

Leaf senescence is one of the programmed cell death (PCD) and is regulated by the ethylene (ET) hormone [35]. Although SA regulates key pathways of plant cell death and immunity, this regulation is inhibited by ET [33]. Interestingly, ET did not affect leaf senescence in young seedlings [36]. There is a sufficient possibility that newly generated upper leaves can escape the effects of ET when used for agro-infiltration. *Arabidopsis* onset of leaf death (*old*) mutants is allelic to *CONSTITUTIVE EXPRESSER OF PR GENES (CPR5)*, which showed senescence symptoms in young seedlings [37, 38]. Interestingly, CPR5 maintains the steady state level of nicotinamide adenine dinucleotide (NAD) [39], suggesting that NAD homeostasis also may affect NLR-mediated cell death including autoimmune

activity. CPR5 associated with a novel nucleoporin PLANT NUCLEAR ENVELOPE TRANSMEMBRANE 1 (PNET1) [40]. Function of human PNET1 homologs is important to cell cycle regulation [41]. This suggests that CPR5-PNET1 may have a dual function between cell cycle and immune pathways in the nuclear pore. Thus, leaf position-dependent HR cell death might be associated with cell cycle and immunity.

Conclusions

Optimized transient HR cell death assay conditions for NLR studies using tobacco plants are proposed. When temperature, humidity, and leaf location conditions are optimally adjusted, various NLR-mediated effectors can induce HR apoptosis and observe an autoimmune phenotype.

Methods

Plant materials and growth conditions

Nicotiana tabacum cv. Petite Gerard and *Nicotiana benthamiana* plants were sown on soil and grown at 22 °C under long day conditions (16 h light/8 h dark) with 55% relative humidity in growth room. For the HR assay, the plants were placed in the growth room for 24 h after infiltration, and then the plants were transferred to the growth chamber to maintain the humidity at 10% or 65%.

Agrobacterium strains and vector constructions

Agrobacterium tumefaciens GV3101 was used in infiltration assays with tobacco leaves. The *Agrobacterium* strain GV3101 competent cells were thawed on ice and added

1 μ L recombinant plant expression vector, then kept on ice for 5 min. Then the mixture was fast frozen in liquid nitrogen for 5 min, followed by an incubation at 37 °C for 5 min. After that, the mixture was kept in ice for 5 min and added 1 mL fresh Luria Broth (LB) liquid medium. After a culture in shaker for 1 h at 28 °C, 200 rpm, 100 μ L cells were plated on an LB agar plate containing rifampicin (25 mg/L) and kanamycin (100 mg/L), and then cultured for 2 days at 28 °C.

Briefly, genomic fragments of full-length RRS1-R and RRS1-S were PCR-amplified from *Arabidopsis thaliana* genomic DNA, accessions Ws-2 and Col-0. The genomic fragments of RRS1 were PCR-amplified with primers containing 4 bp specific overhangs and *Bsa*I recognition sequence and cloned into the pCR8/GW/TOPO (ThermoFisher). The resulting pCR8 constructs were then used for Golden Gate assembly in pICH86988 [42]. C-cCFP, C-nCerulean, and N-cCFP tags were introduced into RRS1 [21]. *35S::RPS4-HA*, *35S::RPS4-Myc*, *35S::RRS1-R-HF*, *35S::AvrRps4-GFP*, *35S::AvrRps4^{E187A}-GFP*, *35S::AvrRps4^{KRVY/AAAA}-GFP*, *35S::TIR^{RPS4}-GFP*, *35S::TIR^{RPS4(SH/AA)}-GFP*, and *35S::TIR^{RRS1(SH/AA)}-GFP* used in this study were described previously [15–17, 21].

Agrobacterium tumefaciens infiltration in tobacco plant

A single positive colony of *Agrobacterium* was inoculated in 3 mL LB liquid medium (25 mg/L rifampicin and 100 mg/L kanamycin) and cultured for 1 day in a shaker at 28 °C, 200 rpm. Then, 100 μ L *Agrobacterium* cells were transferred into 5 mL fresh LB liquid medium supplemented with the above-mentioned antibiotics and cultured overnight at 28 °C, 200 rpm. The *agrobacterium* cells were pelleted by centrifugation (3000 rpm, 10 min) and resuspended with 1 mL infiltration buffer (10 mM MgCl₂, 10 mM MES, pH 5.6). The re-suspended agro-cells were then diluted to OD₆₀₀=0.5 with infiltration buffer. Fully expanded leaves of 4–5 weeks old tobacco plant were syringe-infiltrated in the abaxial region. The leaf-infiltrated plants were dried with paper, and the humidity in the greenhouse was maintained at 55%. After 1 day, the plants were transferred to growth chamber with a 16-h light/8-h dark photoperiod at 22 °C, 65% humidity or transferred to growth chamber at 22 °C, 10% humidity. It is possible to test in a laboratory with a relative humidity of 10% and at 20–22 °C. The HR phenotype was confirmed at 3–5 dpi.

Acknowledgements

The author would like to thank Dr. Panagiotis F. Sarris for Golden Gate cloning the *35S::cCFP-RRS1-S-nCerulean*, *35S::cCFP-RRS1-R-nCerulean*, and *35S::RRS1-R-cCFP*.

Author contributions

SUH designed the research, performed the experiments, analyzed the data, wrote, and approved the final manuscript. The author read and approved the final manuscript

Funding

This project was funded by the Next-Generation BioGreen 21 Program (Project No. PJ01365301) by Rural Development Administration and National Research Foundation of Korea (NRF) grant funded by the Korea government (MEST) (No.18A18297611), Republic of Korea.

Availability of data and materials

All data generated or analysis during this study is included in published article. The materials are available from the corresponding author on reasonable request.

Declarations

Ethics approval and consent to participate

Not applicable.

Consent for publication

Not applicable.

Competing interests

The authors declare that they have no competing interests.

Received: 29 January 2022 Accepted: 21 April 2022

Published: 2 May 2022

References

- Chaloner TM, Gurr SJ, Bebbler DP. Plant pathogen infection risk tracks global crop yields under climate change. *Nat Clim Chang*. 2021;11(8):710–5.
- Velasquez AC, Castroverde CDM, He SY. Plant-pathogen warfare under changing climate conditions. *Curr Biol*. 2018;28(10):R619–34.
- McCann HC, Guttman DS. Evolution of the type III secretion system and its effectors in plant-microbe interactions. *New Phytol*. 2008;177(1):33–47.
- Pruitt RN, Gust AA, Nurnberger T. Plant immunity unified. *Nat Plants*. 2021;7(4):382–3.
- Bentham A, Burdett H, Anderson PA, Williams SJ, Kobe B. Animal NLRs provide structural insights into plant NLR function. *Ann Bot*. 2017;119(5):827–702.
- Adachi H, Derevnina L, Kamoun S. NLR singletons, pairs, and networks: evolution, assembly, and regulation of the intracellular immunoreceptor circuitry of plants. *Curr Opin Plant Biol*. 2019;50:121–31.
- Narusaka M, Shirasu K, Noutoshi Y, Kubo Y, Shiraishi T, Iwabuchi M, Narusaka Y. RRS1 and RPS4 provide a dual Resistance-gene system against fungal and bacterial pathogens. *Plant J*. 2009;60(2):218–26.
- Birker D, Heidrich K, Takahara H, Narusaka M, Deslandes L, Narusaka Y, Reymond M, Parker JE, O'Connell R. A locus conferring resistance to *Colletotrichum higginsianum* is shared by four geographically distinct *Arabidopsis* accessions. *Plant J*. 2009;60(4):602–13.
- Guo H, Ahn HK, Sklenar J, Huang J, Ma Y, Ding P, Menke FLH, Jones JDG. Phosphorylation-Regulated activation of the *Arabidopsis* RRS1-R/RPS4 immune receptor complex reveals two distinct effector recognition mechanisms. *Cell Host Microbe*. 2020;27(5):769–781 e766.
- Menna A, Nguyen D, Guttman DS, Desveaux D. Elevated temperature differentially influences effector-triggered immunity outputs in *Arabidopsis*. *Front Plant Sci*. 2015;6:995.
- Kiraly L, Hafez YM, Fodor J, Kiraly Z. Suppression of tobacco mosaic virus-induced hypersensitive-type necrotization in tobacco at high temperature is associated with downregulation of NADPH oxidase and superoxide and stimulation of dehydroascorbate reductase. *J Gen Virol*. 2008;89(Pt 3):799–808.
- Desaint H, Aoun N, Deslandes L, Vaillau F, Roux F, Berthome R. Fight hard or die trying: when plants face pathogens under heat stress. *New Phytol*. 2021;229(2):712–34.
- Cesari S, Kanzaki H, Fujiwara T, Bernoux M, Chalvon V, Kawano Y, Shimamoto K, Dodds P, Terauchi R, Kroj T. The NB-LRR proteins RGA4 and RGA5 interact functionally and physically to confer disease resistance. *EMBO J*. 2014;33(17):1941–59.

14. Maqbool A, Saitoh H, Franceschetti M, Stevenson CE, Uemura A, Kanzaki H, Kamoun S, Terauchi R, Banfield MJ. Structural basis of pathogen recognition by an integrated HMA domain in a plant NLR immune receptor. *Elife*. 2015;4:e08709.
15. Huh SU, Cevik V, Ding P, Duxbury Z, Ma Y, Tomlinson L, Sarris PF, Jones JDG. Protein-protein interactions in the RPS4/RRS1 immune receptor complex. *PLoS Pathog*. 2017;13(5): e1006376.
16. Williams SJ, Sohn KH, Wan L, Bernoux M, Sarris PF, Segonzac C, Ve T, Ma Y, Saucet SB, Ericsson DJ, et al. Structural basis for assembly and function of a heterodimeric plant immune receptor. *Science*. 2014;344(6181):299–303.
17. Sohn KH, Hughes RK, Piquerez SJ, Jones JD, Banfield MJ. Distinct regions of the *Pseudomonas syringae* coiled-coil effector AvrRps4 are required for activation of immunity. *Proc Natl Acad Sci USA*. 2012;109(40):16371–6.
18. Tasset C, Bernoux M, Jauneau A, Pouzet C, Briere C, Kieffer-Jacquino S, Rivas S, Marco Y, Deslandes L. Autoacetylation of the *Ralstonia solanacearum* effector PopP2 targets a lysine residue essential for RRS1-R-mediated immunity in *Arabidopsis*. *PLoS Pathog*. 2010;6(11): e1001202.
19. Grund E, Tremousaygue D, Deslandes L. Plant NLRs with integrated domains: unity makes strength. *Plant Physiol*. 2019;179(4):1227–35.
20. Ma Y, Guo H, Hu L, Martinez PP, Moschou PN, Cevik V, Ding P, Duxbury Z, Sarris PF, Jones JDG. Distinct modes of derepression of an *Arabidopsis* immune receptor complex by two different bacterial effectors. *Proc Natl Acad Sci USA*. 2018;115(41):10218–27.
21. Sarris PF, Duxbury Z, Huh SU, Ma Y, Segonzac C, Sklenar J, Derbyshire P, Cevik V, Rallapalli G, Saucet SB, et al. A plant immune receptor detects pathogen effectors that target WRKY transcription factors. *Cell*. 2015;161(5):1089–100.
22. Maruta N, Burdett H, Lim BYJ, Hu X, Desa S, Manik MK, Kobe B. Structural basis of NLR activation and innate immune signalling in plants. *Immunogenetics*. 2022;74:5–26.
23. Essuman K, Summers DW, Sasaki Y, Mao X, Yim AKY, DiAntonio A, Milbrandt J. TIR domain proteins are an ancient family of NAD(+)-Consuming enzymes. *Curr Biol*. 2018;28(3):421–430 e424.
24. Horsefield S, Burdett H, Zhang X, Manik MK, Shi Y, Chen J, Qi T, Gilley J, Lai JS, Rank MX, et al. NAD(+) cleavage activity by animal and plant TIR domains in cell death pathways. *Science*. 2019;365(6455):793–9.
25. Wan L, Essuman K, Anderson RG, Sasaki Y, Monteiro F, Chung EH, Osborne Nishimura E, DiAntonio A, Milbrandt J, Dangl JL, et al. TIR domains of plant immune receptors are NAD(+)-cleaving enzymes that promote cell death. *Science*. 2019;365(6455):799–803.
26. Li Q, Jiang XM, Shao ZQ. Genome-wide analysis of NLR disease resistance genes in an updated reference genome of Barley. *Front Genet*. 2021;12: 694682.
27. Wang L, Zhao L, Zhang X, Zhang Q, Jia Y, Wang G, Li S, Tian D, Li WH, Yang S. Large-scale identification and functional analysis of NLR genes in blast resistance in the Tetep rice genome sequence. *Proc Natl Acad Sci U S A*. 2019;116(37):18479–87.
28. Van de Weyer AL, Monteiro F, Furzer OJ, Nishimura MT, Cevik V, Witek K, Jones JDG, Dangl JL, Weigel D, Bemm F. A Species-Wide inventory of NLR Genes and Alleles in *Arabidopsis thaliana*. *Cell*. 2019;178(5):1260–1272 e1214.
29. Huang Z, Qiao F, Yang B, Liu J, Liu Y, Wulff BBH, Hu P, Lv Z, Zhang R, Chen P, et al. Genome-wide identification of the NLR gene family in *Haynaldia villosa* by SMRT-RenSeq. *BMC Genomics*. 2022;23(1):118.
30. Stam R, Scheikl D, Tellier A. Pooled enrichment sequencing identifies diversity and evolutionary pressures at NLR Resistance Genes within a wild tomato population. *Genome Biol Evol*. 2016;8(5):1501–15.
31. Jupe F, Chen X, Verweij W, Witek K, Jones JD, Hein I. Genomic DNA library preparation for resistance gene enrichment and sequencing (RenSeq) in plants. *Methods Mol Biol*. 2014;1127:291–303.
32. Disch EM, Tong M, Kotur T, Koch G, Wolf CA, Li X, Hoth S. Membrane-Associated Ubiquitin Ligase SAUL1 suppresses temperature- and humidity-dependent autoimmunity in *Arabidopsis*. *Mol Plant Microbe Interact*. 2016;29(1):69–80.
33. Li Z, Liu H, Ding Z, Yan J, Yu H, Pan R, Hu J, Guan Y, Hua J. Low temperature enhances plant immunity via salicylic acid pathway genes that are repressed by Ethylene. *Plant Physiol*. 2020;182(1):626–39.
34. Zhang Y, Goritschnig S, Dong X, Li X. A gain-of-function mutation in a plant disease resistance gene leads to constitutive activation of downstream signal transduction pathways in suppressor of npr1-1, constitutive 1. *Plant Cell*. 2003;15(11):2636–46.
35. Koyama T. A hidden link between leaf development and senescence. *Plant Sci*. 2018;276:105–10.
36. Jing HC, Schippers JH, Hille J, Dijkwel PP. Ethylene-induced leaf senescence depends on age-related changes and OLD genes in *Arabidopsis*. *J Exp Bot*. 2005;56(421):2915–23.
37. Jing HC, Dijkwel PP. CPR5: a Jack of all trades in plants. *Plant Signal Behav*. 2008;3(8):562–3.
38. Jing HC, Hebel R, Oeljeklaus S, Sitek B, Stuhler K, Meyer HE, Sturre MJ, Hille J, Warscheid B, Dijkwel PP. Early leaf senescence is associated with an altered cellular redox balance in *Arabidopsis* cpr5/old1 mutants. *Plant Biol (Stuttg)*. 2008;10(Suppl 1):85–98.
39. Schippers JH, Nunes-Nesi A, Apetrei R, Hille J, Fernie AR, Dijkwel PP. The *Arabidopsis* onset of leaf death5 mutation of quinolinate synthase affects nicotinamide adenine dinucleotide biosynthesis and causes early ageing. *Plant Cell*. 2008;20(10):2909–25.
40. Tang Y, Huang A, Gu Y. Global profiling of plant nuclear membrane proteome in *Arabidopsis*. *Nat Plants*. 2020;6(7):838–47.
41. Fujitomo T, Daigo Y, Matsuda K, Ueda K, Nakamura Y. Critical function for nuclear envelope protein TMEM209 in human pulmonary carcinogenesis. *Cancer Res*. 2012;72(16):4110–8.
42. Weber E, Engler C, Gruetzner R, Werner S, Marillonnet S. A modular cloning system for standardized assembly of multigene constructs. *PLoS ONE*. 2011;6(2): e16765.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

